

# NeoNatalie

Realistic and affordable training and therapy



In support of the UN Millennium Development Goal 4 Laerdal has committed to providing the NeoNatalie Newborn Simulator and the NeoNatalie resuscitation tools to developing countries on a not-for-profit basis through 2015





**NeoNatalie** facilitates effective learning of standard newborn care and resuscitation measures. The low cost simulator is available with affordable therapeutic tools.



NeoNatalie is easily prepared for use by filling the body with two liters of water, or alternatively with air.



NeoNatalie facilitates practice in effective bag-mask-ventilation. Chest will only rise with correct technique.



NeoNatalie was developed to help train millions of birth attendants in low-resource countries in neonatal resuscitation courses, such as “**Helping Babies Breathe**”\*

\* Helping Babies Breathe (HBB) is an initiative of the American Academy of Pediatrics in collaboration with the World Health Organization (WHO), US Agency for International Development, Saving Newborn Lives, the National Institute of Child Health and Human Development (NICHD), and a number of other global health organizations.

For more information, see [www.helpingbabiesbreathe.org](http://www.helpingbabiesbreathe.org)



## Highly realistic simulation training at a low cost

NeoNatalie has realistic size and appearance, and also natural weight, feel and touch when filled with lukewarm water. The baby’s status can be simulated as desired for effective role playing in relevant scenarios covering normal care after birth as well as standard resuscitation measures.



Use squeeze bulbs to simulate

- birth cries
- spontaneous breathing
- palpable umbilical pulse
- heart sounds



## NeoNatalie Resuscitator

The bag-valve-mask unit is made of silicone and polysulphonate, can be boiled or autoclaved, and is extremely durable. Comes with two mask sizes. Time proven design and quality. An oxygen kit is an optional extra.

## NeoNatalie Suction

The suction unit is made in one piece of silicone, can be boiled or autoclaved, and withstands hundreds of uses. Shaped like a penguin, its beak is ideal for newborn oral and nasal suction. Head can be tilted aside for easy emptying and cleaning.

## ORDER INFORMATION

NeoNatalie Newborn Simulator is available with a light or dark complexion in two models.

### Basic Model

with Standard Accessories

- Simulation squeeze bulbs with tubing and connectors
- Umbilical cord with connector and two ties
- Two sheets to simulate towels
- Head cap
- Tube for topping body filling
- Directions for use
- Transport/storage bag



### Complete Model

with additional resuscitation care kit

- NeoNatalie Resuscitator
- NeoNatalie Suction
- Training stethoscope



NeoNatalie comes with a shoulder bag for convenient transport and storage

Cat. no.	Title	Includes
104-10001	NeoNatalie Basic (dark)	NeoNatalie Newborn Simulator; dark complexion w/ standard accessories
104-10002	NeoNatalie Basic (light)	NeoNatalie Newborn Simulator; light complexion w/ standard accessories
104-10005	NeoNatalie Complete (dark)	NeoNatalie Newborn Simulator; dark complexion w/ standard accessories and resuscitation care kit
104-10006	NeoNatalie Complete (light)	NeoNatalie Newborn Simulator; light complexion w/ standard accessories and resuscitation care kit
846040	NeoNatalie Resuscitator	Resuscitator; 1 mask size 0, 1 mask size 1, directions for use
986000	NeoNatalie Suction	Bulb suction, directions for use

To order and for more information, see [www.laerdalglobalhealth.com](http://www.laerdalglobalhealth.com)