



1. Allow the cup to air dry in a protected space.
2. Inspect for damage and cleanliness/mineral deposits.
3. If damaged, remove cup from use.

STEP 3




1. Use a fresh solution of 0.5% chlorine.
2. Immerse cup completely in the solution for 20 min.
3. Shake off and empty chlorine solution from the cup.
4. Rinse the cup thoroughly in clean water.
5. Repeat 3 times with fresh water to remove any chemical residue.

c) Chlorine (0.5%)



136 °C (276 °F)
2.0 kg/cm²
10-20 minutes

b) Steam Autoclaving




100 °C (212 °F)
10 minutes.

a) Boiling

Perform one of the listed methods:

STEP 2



1. Wash in a clean tray with clean water and mild soap.
2. Scrub the cup to remove milk and all other residue.
3. Rinse in clean water to remove soap residue.
4. Repeat above steps until clean.

STEP 1

Reprocessing Instructions

EN

Intended Use

Nifty Feeding Cup is a reusable cup for feeding breast milk to newborns with breastfeeding difficulties.

- ⚠ Caution:**
- Reprocess the cup before first use and with each new infant. Refer to Reprocessing Instructions.
 - Do not use the cup to administer medication.


Use

Express desired amount of breast milk in a clean Nifty cup.

Gently squeeze the cup to fill the feeding reservoir.

Ensure baby is in a semi-upright position with a supported head and neck.

Bring the cup to baby's mouth. Do not pour the milk into the mouth. Allow baby to actively feed from the reservoir at her own pace.



Design Collaboration

Nifty Feeding Cup has been developed and designed through a collaboration between the University of Washington, Seattle Children's Hospital, PATH and Laerdal Global Health.

© 2016 Laerdal Medical AS, All rights reserved.

Laerdal® is a registered trademark of Laerdal Medical AS.

Manufactured in China for:

Laerdal Medical AS
PO. Box 377 Tanke Svlansgate 30, 4002 Stavanger, Norway
Tel : +47 51 51 17 00
Fax: +47 51 52 35 57
Email: LGH@laerdal.com

Printed in China

Nifty Feeding Cup

