



Moyo

## Fetal heart rate monitor

A compact digital fetal heart rate monitor used during labor to help health workers make timely decisions by providing accurate feedback on fetal heart rate.



# Moyo

## Accurate fetal heart rate monitoring

### Timely and accurate monitoring

A 9-crystal sensor accurately detects the fetal heart rate as well as differentiates the maternal from the fetal heart rate.

### Aids fast reaction to abnormal fetal heart rate

Early detection of abnormal FHR guides decision making and helps health workers to act quickly. When used continuously, an audio-visual alarm is activated if an abnormal FHR is detected, and a 30 min histogram is provided.

### Continuous and intermittent monitoring

Moyo is lightweight and portable – it can monitor continuously, even while the mother is moving around. It eases the workload for health workers without interrupting existing routines.

## Features

- Both intermittent and continuous monitoring
- 9-crystal sensor for better accuracy
- Differentiates the maternal and the fetal heart rate
- 30-minute histogram display of the FHR
- Audio-visual alarm if abnormal FHR is detected
- Lightweight and portable, can monitor continuously, even while the mother is moving around
- Internal rechargeable battery



The sensor can be strapped on to the mother's abdomen. Lightweight and portable, Moyo allows the mother to move around freely.



Normal fetal heart rate



Abnormal fetal heart rate



Prolonged Abnormal fetal heart rate



30- minute histogram