



NeoNatalie Live

Smart newborn resuscitation trainer

Built to promote scalable and high-quality training, NeoNatalie Live empowers healthcare providers and students to practice newborn resuscitation skills, frequently and independently.



NeoNatalie Live

Smart Newborn Resuscitation Trainer

Autonomous, high-quality training

NeoNatalie Live can act as a skills corner where learners engage in realistic scenarios based on 1,200 real ventilation cases at their own pace, supporting frequent, high-quality practice that builds confidence and competence.

Remote monitoring and scalable feedback

Training data is collected and accessible remotely, allowing supervisors to monitor progress and identify areas for improvement. This makes it possible to scale training across locations while maintaining consistent quality and oversight.

Consistent training across all care levels

From tertiary hospitals to primary health clinics, as well as educational institutions, NeoNatalie Live supports WHO's Essential Newborn Care Course and enables equitable access to trainings in newborn resuscitation.

Plays a pivotal role for service providers

"Where time and training opportunities are limited for healthcare personnel, the NeoNatalie Live enables self-directed sessions for individuals or teams with the ongoing objective assessment and feedback. NeoNatalie Live is playing a pivotal role for the service providers in development of skill, efficiency and promoting self-evaluation and knowledge retention in case of newborn resuscitation and survival"

Tamanaa Eshika,
Project Manager for SWAP

Features

- Short and self-directed sessions
- Four different patient cases
- Based on over 1,200 resuscitations
- Objective feedback on performance
- Detection of head tilt
- Visible chest rise
- Ventilation performance impacts heart rate
- Crying sounds when successful
- Compatible with NeoBeat newborn heart rate meter

